

## 4-Primer Amplicon Tagging with the LP 48.48 IFC


### IMPORTANT

- Before using this document, read and understand the detailed instructions and safety guidelines in the Access Array™ System for Illumina® Sequencing Systems User Guide (PN 100-3770) or the Access Array System for the Ion Torrent™ PGM™ Sequencing System User Guide (PN 100-5024).
- This quick reference requires the use of updated scripts for FC1 (v1.6 or later) and IFC Controller AX (v2.8 or later) that can be downloaded from [fluidigm.com/software](http://fluidigm.com/software).

### Prime the IFC

#### IMPORTANT

- Use the LP 48.48 integrated fluidic circuit (IFC) or the Access Array 48.48 IFC within 24 hours of opening the package.
- Due to different accumulator volumes, only use 48.48 syringes with 300  $\mu\text{L}$  of Control Line Fluid (LP 48.48: PN 100-2345; Access Array 48.48: PN 89000020).
- Control Line Fluid on IFC or in the inlets makes the IFC unusable.
- Load the IFC within 60 minutes of priming.
- Be certain that the 1X Access Array Harvest Solution and 1X Access Array Hydration Reagent v2 are thawed completely to room temperature and mixed thoroughly prior to use.



**NOTE** The location of the sample inlets is different from 48.48 Gene Expression and Genotyping IFCs.

- 1 Inject Control Line Fluid into each accumulator on the IFC.
- 2 Add 500  $\mu\text{L}$  of 1X Access Array Harvest Solution (Fluidigm, PN 100-1031) into the H1–H3 wells of the IFC.
- 3 Add 500  $\mu\text{L}$  of Access Array Hydration Reagent v2 (● blue cap, Fluidigm, PN 100-7966) into the H4 well of the IFC.  
**NOTE** 1X Access Array Harvest Solution (Fluidigm, PN 100-1031) can be substituted for Access Array Hydration Reagent v2 (blue cap, Fluidigm, PN 100-7966) in the H4 well.
- 4 Remove and discard the protective film from bottom of the IFC.
- 5 In the pre-PCR lab, load the IFC into the **pre-PCR** IFC Controller AX.
- 6 Prime the LP 48.48 IFC. Select **Prime (155x)** and **Run Script**.

**NOTE** For Access Array 48.48 IFCs, select **Prime (151x)** and **Run Script**. If your instrument does not contain the Prime (155x) script, download the latest version of the IFC Controller AX instrument software from [fluidigm.com/software](http://fluidigm.com/software).

### Prepare the 20X Primer Solution

**IMPORTANT** Be certain that the 20X Access Array Loading Reagent is thawed completely to room temperature and mixed thoroughly prior to use.

| Component  | Volume ( $\mu\text{L}$ ) |
|--|--------------------------|
| 50 $\mu\text{M}$ CS1-tagged TS forward primer*   | 2.0                      |
| 50 $\mu\text{M}$ CS2-tagged TS reverse primer    | 2.0                      |
| 20X Access Array Loading Reagent (PN 100-0883) ○ | 5.0                      |
| PCR-certified water                              | 91.0                     |
| <b>Total</b>                                     | <b>100.0</b>             |

\*Fluidigm Access Array target-specific primers (PN ASY-AA) are provided as mixed forward and reverse primers at 50  $\mu\text{M}$  per primer. If using Fluidigm assays, pipet 2  $\mu\text{L}$  of mixed primers and adjust water to 93  $\mu\text{L}$ .

Vortex the 20X primer solution for 20 seconds and centrifuge for 30 seconds.

**NOTE** The final concentration of each primer is 1  $\mu\text{M}$  in the 20X primer solution and 50 nM in the PCR reaction.

### Prepare the Sample Pre-Mix Solution

- 1 In a DNA-free hood, combine the following components from the FastStart™ High Fidelity PCR System, dNTPack (Sigma-Aldrich, PN 04738292001), with 20X Access Array Loading Reagent and PCR-certified water in a 1.5 mL sterile tube (sufficient volume for one IFC):

| Component   | Volume ( $\mu\text{L}$ ) | Volume per IFC* ( $\mu\text{L}$ ) |
|---|--------------------------|-----------------------------------|
| PCR-certified water   | 0.95                     | 57.0                              |
| 10X FastStart High Fidelity Reaction Buffer without $\text{MgCl}_2$ | 0.5                      | 30.0                              |
| 25 mM $\text{MgCl}_2$   | 0.9                      | 54.0                              |
| DMSO  | 0.25                     | 15.0                              |
| 10 mM PCR Grade Nucleotide Mix                                      | 0.1                      | 6.0                               |
| 20X Access Array Loading Reagent (PN 100-0883) ○                    | 0.25                     | 15.0                              |
| 5 U/ $\mu\text{L}$ FastStart High Fidelity Enzyme Blend             | 0.05                     | 3.0                               |
| <b>Total</b>  | <b>3.0</b>               | <b>180.0</b>                      |

\*60 reactions for ease of pipetting into a single LP 48.48 IFC or Access Array 48.48 IFC

- 2 Vortex the sample pre-mix solution for 20 seconds and centrifuge for 30 seconds before preparing the sample mix solution.

## Prepare the Sample Mix Solution

Primers containing a unique, sample-specific barcode must be added to each sample.

- Combine the following components in a 96-well PCR plate to prepare 48 individual sample mix solutions where each well receives a uniquely barcoded primer pair:

| Component                          | Volume (μL) |
|------------------------------------|-------------|
| Sample pre-mix solution            | 3.0         |
| 50 ng/μL genomic DNA               | 1.0         |
| 2 μM Access Array Barcode Primers* | 1.0         |
| <b>Total</b>                       | <b>5.0</b>  |

\*Select 48 unique barcodes from either library: Access Array Barcode Library for Illumina Sequencers—384 (Single direction; PN 100-4876) or Access Array Barcode Library for the Ion Torrent PGM™ Sequencer—96 (Bidirectional; PN 100-4911).


- Vortex for 20 seconds and centrifuge for 30 seconds.

## Load the IFC


### IMPORTANT

- Check IFC orientation before pipetting reagents into inlets.
- While pipetting, do not go past the first stop on the pipette. Doing so might introduce air bubbles into the inlets.

### Loading Map for LP 48.48 IFC or Access Array 48.48 IFC



### Pipetting Scheme



- Pipet 4 μL of 20X primer solution into each primer inlet.
- Pipet 4 μL of sample mix solution into each sample inlet.

**IMPORTANT** All inlets must be filled. If you have less than 48 primer sets or samples, prepare mixes as follows: Use water in place of tagged primers in the 20X primer solution and use water in place of the genomic DNA and barcode primers in the sample mix solution.

- Load the LP 48.48 IFC. Select **Load Mix (155x)** and **Run Script**.

**NOTE** For Access Array 48.48 IFCs, select **Load Mix v7 (151x)** and **Run Script**. Load Mix v7 (151x) is a script update from Load Mix (151x). If your instrument does not contain the Load Mix v7 (151x) script, download the latest version of the IFC Controller AX instrument software from [fluidigm.com/software](http://fluidigm.com/software).

## Thermal-Cycle the IFC


Choose the thermal-cycling instrument and run the protocol:

| Instrument                         | Run Protocol                |
|------------------------------------|-----------------------------|
| FC1™ cyclers*                      | <b>AA 48x48 Standard v1</b> |
| Stand-alone thermal cycler (SATC)† | <b>AA48v1</b>               |

\*See the FC1 Cycler Usage Quick Reference (100-1250).

†See the Stand-Alone Thermal Cycler Quick Reference (68000111).

## Harvest the IFC



- In the post-PCR lab, replace H1–H4 with 650 μL of 1X Access Array Harvest Solution. (Do not use the Hydration Reagent v2.)
- Pipet 2 μL of 1X Access Array Harvest Solution into each of the sample inlets on the IFC.
- Load the IFC into the **post-PCR** IFC Controller AX.
- Harvest the LP 48.48 IFC. Select **Harvest (155x)** and **Run Script**.

### NOTE

• For Access Array 48.48 IFCs, select **Harvest v7 (151x)** and **Run Script**. Harvest v7 (151x) is a script update from Harvest v5 (151x). If your instrument does not contain the Harvest v7 (151x) script, download the latest version of the IFC Controller AX instrument software from [fluidigm.com/software](http://fluidigm.com/software).

• Label a 96-well plate using the barcode number on the IFC. Transfer the harvested samples into columns 1–6 of the pre-labeled 96-well PCR plate. See **Pipetting Scheme**.

- Continue with your library preparation as described in the Post-PCR Amplicon Purification and Quantitation chapter of the Illumina (PN 100-3770) or Ion Torrent (PN 100-5024) user guide.

## For technical support visit [fluidigm.com/support](http://fluidigm.com/support).

North America +1 650 266 6100 | Toll-free (US/CAN): 866 358 4354 | [techsupport@fluidigm.com](mailto:techsupport@fluidigm.com)  
 Europe/Middle East/Africa/Russia +44 1223 598100 | [techsupport@fluidigm.com](mailto:techsupport@fluidigm.com)  
 Japan +81 3 3662 2150 | [techsupport@fluidigm.com](mailto:techsupport@fluidigm.com)

Central and South America +1 650 266 6100 | [techsupport@fluidigm.com](mailto:techsupport@fluidigm.com)  
 China (excluding Hong Kong) +86 21 3255 8368 | [techsupportchina@fluidigm.com](mailto:techsupportchina@fluidigm.com)  
 All other Asian countries/India/Australia +1 650 266 6100 | [techsupportasia@fluidigm.com](mailto:techsupportasia@fluidigm.com)

## For Research Use Only. Not for use in diagnostic procedures.

Information in this publication is subject to change without notice. Safety data sheet information: [fluidigm.com/sds](http://fluidigm.com/sds). Patent and license information: [fluidigm.com/legalnotices](http://fluidigm.com/legalnotices). EU's WEEE directive information: [fluidigm.com/compliance](http://fluidigm.com/compliance). Fluidigm, the Fluidigm logo, Access Array, and FC1 are trademarks and/or registered trademarks of Fluidigm Corporation in the United States and/or other countries. All other trademarks are the sole property of their respective owners. © 2018 Fluidigm Corporation. All rights reserved. 08/2018