

C1 Single-Cell mRNA Seq HT Reagent Kit v2

Kit Part Number: 101-3473








Purpose

This kit provides the necessary materials to enable sample preparation for mRNA sequencing with high-throughput integrated fluidic circuits (HT IFCs) on the C1™ system.

Contents

When ordering the modules from Fluidigm, use the kit part number: 101-3473.

Part Number	Component	Cap Color	Unit Quantity	Volume/Unit	Uses/Unit	Storage
101-3474	C1 Single-Cell mRNA Seq HT Reagent Kit v2, Module 1					
101-0165	10X Blocking Reagent		5 tubes	50.0 µL	1 IFC	4 °C
100-8506	Cell Rinsing Reagent		7 tubes	1.6 mL	1 IFC	
101-0552	Stability Solution		1 tube	1.0 mL	5 IFCs	
100-5315	Suspension Reagent		5 tubes	50.0 µL	1 IFC	
101-2132	Valve Fluid V2		3 tubes	1.5 mL	1 IFC	
101-3334	Actuation Fluid		2 tubes	1.5 mL	1 IFC	
101-0060	C1 Single-Cell mRNA Seq HT Reagent Kit, Module 2					
100-6248	C1 Harvest Reagent		5 bottles	3.7 mL	1 IFC	-20 °C
100-5361	DNA Dilution Reagent		1 tube	1.0 mL	5 IFCs	
101-0058	Enrichment Primer		1 tube	130.0 µL	5 IFCs	
100-5312	Loading Reagent		5 tubes	25.0 µL	1 IFC	
101-0054	Preamp Primer		1 tube	15.0 µL	5 IFCs	
101-1125	Preloading Reagent		5 tubes	750.0 µL	1 IFC	
101-0061	C1 Single-Cell mRNA Seq HT Reagent Kit, Module 3					
101-0164	Cell Barcode Plate	—	5 plates	40 barcodes	Single use only	-20 °C
101-0062	C1 Single-Cell mRNA Seq HT Reagent Kit, Module 4					
100-0056	Reverse Transcription Primer		1 tube	20.0 µL	5 IFCs	-80 °C

Part Number	Component	Cap Color	Unit Quantity	Volume/Unit	Uses/Unit	Storage
101-5823	C1 Single-Cell mRNA Seq HT Reagent Kit v2, Module 5					
	5X First-Strand Buffer		1 tube	96.0 µL	5 IFCs	-20 °C
	dNTP Mix (20 mM each)		1 tube	30.0 µL	5 IFCs	
	Dithiothreitol (DTT; 100mM)		1 tube	15.0 µL	5 IFCs	
	Reverse Transcriptase (100U/µL)		1 tube	48.0 µL	5 IFCs	
	Nuclease-Free Water		1 tube	1.0 mL	5 IFCs	
	RNase Inhibitor (40U/µL)		1 tube	60.0 µL	5 IFCs	
	10X Lysis Buffer – v3		1 tube	260.0 µL	5 IFCs	

Handling

Before using reagents, refer to the Generate cDNA Libraries with the C1 Single-Cell mRNA Seq HT IFC and Reagent Kit v2 Protocol (PN 101-4964).

Safety

Use standard laboratory safety protocols. Read and understand the safety data sheets (SDSs) before handling chemicals. To obtain SDSs, go to fluidigm.com/sds and search for the SDS using either the product name or the part number.

For technical support visit fluidigm.com/support.

North America +1 650 266 6100 | Toll-free (US/CAN): 866 358 4354 | techsupport@fluidigm.com **Latin America** +1 650 266 6100 | techsupportlatam@fluidigm.com
Europe/Middle East/Africa/Russia +44 1223 859941 | techsupporteurope@fluidigm.com **China (excluding Hong Kong)** +86 21 3255 8368 | techsupportchina@fluidigm.com
Japan +81 3 3662 2150 | techsupportjapan@fluidigm.com **All other Asian countries/India/Australia** +1 650 266 6100 | techsupportasia@fluidigm.com

For Research Use Only. Not for use in diagnostic procedures.

Information in this publication is subject to change without notice. **Safety data sheet information:** fluidigm.com/sds. **Patent and license information:** fluidigm.com/legalnotices. **EU's WEEE directive information:** fluidigm.com/compliance. Fluidigm, the Fluidigm logo, and C1 are trademarks or registered trademarks of Fluidigm Corporation in the United States and/or other countries. All other trademarks are the sole property of their respective owners.
 © 2017 Fluidigm Corporation. All rights reserved. 06/2017