

# Cell-ID Intercalator-Ir—500 $\mu\text{M}$

Catalog: 201192B  
Package size: 500  $\mu\text{L}$

## Storage:

- Upon receiving this product, divide it into aliquots and freeze them at  $-20\text{ }^{\circ}\text{C}$ .
- Aliquots stored at  $4\text{ }^{\circ}\text{C}$  are stable for up to three months.
- Frozen aliquots should be used only once after thawing.



**WARNING** Before handling any chemicals, refer to the safety data sheet (SDS) provided by the manufacturer, and observe all relevant precautions.

## Description

Cell-ID™ Intercalator-Ir is a cationic nucleic acid intercalator that contains natural abundance iridium ( $^{191}\text{Ir}$  and  $^{193}\text{Ir}$ ) and is used for identifying nucleated cells in CyTOF® system analysis. When cells are stained with Intercalator-Ir, it binds to cellular nucleic acid, and detection of both stable isotopes enables identification of nucleated cells. It is a live-cell membrane-impermeable dye and therefore requires cells to be fixed and/or permeabilized before staining.

## Important Product Notes

- Cell-ID Intercalator-Ir 500  $\mu\text{M}$  is a highly concentrated metal intercalator solution. It must be diluted in accordance with this protocol to avoid early failure of the detector.
- While dilutions of the 500  $\mu\text{M}$  stock solution are suggested in the protocol below, the concentration can be titrated for individual cell types and experiments for optimal Cell-ID Intercalator staining. It is suggested that the intercalator concentration in the staining solution not exceed 1  $\mu\text{M}$ .

## Staining Protocol

- 1 After cell staining is complete, prepare 1 mL of cell intercalation solution for each sample by diluting Cell-ID Intercalator-Ir 1:4,000 into Maxpar® Fix and Perm Buffer (Fluidigm Cat. 201067) and mix by vortexing.
- 2 Add 1 mL of the intercalation solution prepared in step 1 to each tube and gently vortex. Incubate for one hour at room temperature or leave overnight at  $4\text{ }^{\circ}\text{C}$ .

**Note:** Cells can be left at  $4\text{ }^{\circ}\text{C}$  in the intercalation solution up to 48 hours.

- 3 Wash cells by adding 2 mL of Maxpar Cell Staining Buffer (Cat. 201068), centrifuge and discard supernatant by aspiration.
- 4 Repeat for a total of two washes with Maxpar Cell Staining Buffer.
- 5 Wash cells with 2 mL of Maxpar Water (Cat. 201069), centrifuge and discard supernatant by aspiration.
- 6 Leave cells pelleted until ready to run on the CyTOF system. Immediately prior to data acquisition, adjust cell concentration to  $2.5\text{--}5 \times 10^5/\text{mL}$  with Maxpar Water and filter cells into cell strainer cap tubes.
- 7 Acquire data on the CyTOF system.

**For technical support visit [www.fluidigm.com/support](http://www.fluidigm.com/support).**

**North America** +1 650 266 6100 | Toll-free (US/CAN): 866 358 4354 | [support.northamerica@fluidigm.com](mailto:support.northamerica@fluidigm.com) **Latin America** +1 650 266 6100 | [techsupportlatam@fluidigm.com](mailto:techsupportlatam@fluidigm.com)  
**Europe/Middle East/Africa/Russia** +44 1223 859941 | [support.europe@fluidigm.com](mailto:support.europe@fluidigm.com) **China (excluding Hong Kong)** +86 21 3255 8368 | [techsupportchina@fluidigm.com](mailto:techsupportchina@fluidigm.com)  
**Japan** +81 3 3662 2150 | [techsupportjapan@fluidigm.com](mailto:techsupportjapan@fluidigm.com) **All other Asian countries/India/Australia** +1 650 266 6100 | [techsupportasia@fluidigm.com](mailto:techsupportasia@fluidigm.com)

**For Research Use Only. Not for use in diagnostic procedures.**

Information in this publication is subject to change without notice. **Safety data sheet information:** [fluidigm.com/sds](http://fluidigm.com/sds). **Patent and license information:** [fluidigm.com/legalnotices](http://fluidigm.com/legalnotices). **EU's WEEE directive information:** [fluidigm.com/compliance](http://fluidigm.com/compliance). Fluidigm, the Fluidigm logo, Cell-ID, CyTOF, and Maxpar are trademarks or registered trademarks of Fluidigm Corporation in the United States and/or other countries. © 2017 Fluidigm Corporation. All rights reserved. 01/2017