

# DELTAgene™ Assays

Part Numbers: ASY-GE or ASY-GE-WET

## Intended Use

Gene expression assays are intended for use with BioMark™ HD System, reagents, and IFCs (sold separately). Please refer to the BioMark™ HD System User Guide, PN 100-2451, for protocol references.

Contents				
PART NUMBERS	REAGENT	CONCENTRATION of each individual primer*	VOLUME	STORAGE
ASY-GE or ASY-GE-WET	Mixed Forward and Reverse Primers	100 µM	100 µl	-20°C

\*Assay plate contains 50 µL of 200 µM forward mixed with 50 µL of 200 µM reverse primers.

## Handling:

Before using reagents:

1. Bring to room temperature.
2. Centrifuge for 30 seconds to collect materials.
3. Remove original plate seal. Original seal is intended for single use only.
4. Manually mix the contents of the plate. (Do not vortex to avoid cross contamination)
5. Reseal the plate with a new seal.
6. Centrifuge for 30 seconds to collect contents.

## Safety:

Use standard laboratory safety protocols. For MSDS information, contact Fluidigm Technical Support.

## Technical Support:

United States: 1-866-FLUIDLINE (1.866.358.4354) or techsupport@fluidigm.com

Outside the United States: 650-266-6100 | [www.fluidigm.com/support.htm](http://www.fluidigm.com/support.htm)

Fluidigm, the Fluidigm logo, BioMark, and DELTAgene are trademarks or registered trademarks of Fluidigm Corporation in the U.S. and/or other countries. Fluidigm recommends you only use licensed PCR assay reagents from authorized sources. ©Fluidigm Corporation. All rights reserved. For research use only. Not for diagnostic use.

### NOTICE TO PURCHASER: Disclaimer of License for Pre-Amp Methods

No right to perform Life Technologies Corporation's patented pre-amplification methods is conveyed with the purchase of this DELTAgene Assay from Fluidigm Corporation. A license to use this product with such pre-amplification methods can be obtained (i) with purchase of a BMK-M-96.96 or BMK-M-48.48 Dynamic Array IFC from Fluidigm Corporation or (ii) by a separate license from Life Technologies Corporation. No right to resell this product and no other rights (such as real-time PCR methods, apparatus, reagents or software to perform digital PCR methods) are conveyed by Life Technologies Corporation expressly, by implication, or by estoppel. For information on obtaining additional rights, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.