

Minimum Specifications for Imaging Cells in Fluidigm Integrated Fluidic Circuits

A moderate- or high-resolution microscope system allows imaging of cells in Fluidigm integrated fluidic circuits (IFCs). A microscope is recommended for confirming capture, viability, and morphology of cells in all cell-based IFCs. The following is a list of recommended specifications for appropriate IFC imaging systems. To perform imaging of cells in IFCs, confirm that a microscope in your lab meets these specifications.

Minimum Specifications for Imaging

Need	Recommended Microscope Setup	
Highest-confidence sensing with automated scanning	 Inverted compound microscope with phase contrast, fluorescence capabilities, and motion-control including: 10X phase contrast objective and matching transmitted light condenser phase ring 10X plan apochromat objective 20X plan apochromat objective (for sub-cellular resolution staining) NOTE Any magnification >20X requires long working distance objectives. Automated motorized stage (XYZ) with positional accuracy better than 1.0 μm and an adapter for standard cell culture plates Auto-focus function that stabilizes images along the Z-axis All associated fluorescence imaging hardware: Filter cubes—GFP and Texas Red® minimally recommended High-intensity fluorescence source 1M pixel or greater camera with 5e-dark noise or better and capable of full-frame capture at 10 Hz. Workstation computer with dedicated graphics card (Preferred) 1,920 x 1,200-capable monitor or better Image processing software capable of stitching "200 high resolution images 	
Highest-confidence sensing with manual scanning	 Inverted compound microscope with phase contrast and fluorescence capabilities including: 10X phase contrast objective and matching transmitted light condenser phase ring 10X plan apochromat objective All associated fluorescence imaging hardware: Filter cubes—GFP and Texas Red minimally recommended High-intensity fluorescence source (Preferred) 1M pixel or greater camera with 5e-dark noise or better and capable of full-frame capture at 10 Hz (Preferred) Workstation computer with dedicated graphics, 1,920 x 1,200-capable monitor or better, and image processing software 	
Moderate-confidence sensing with manual scanning	 In order from highest to lowest confidence: Inverted compound microscope with phase contrast capability with 10X phase contrast objective and matching transmitted light condenser phase ring Upright compound microscope with DIC and 10X objective 	

Suggested Fluorescence Microscopes

The table is a list of suggest microscopes that have been tested by the Fluidigm R&D team.

Product Name	Company	Part Number
Automated Inverted Microscope	Leica®	DMi8
Axio Observer Z1*	ZEISS	491912-0003-000

*The optional Perfect Focus module for image stabilization along the Z-axis is recommended for faster imaging.

Fluidigm Solutions for Single-Cell and Cell Biology Applications

Product	Description
C1 [™] Single-Cell Open App [™] IFCs	IFCs for single-cell capture, lysis, and processing on the C1 system
C1 Single-Cell Preamp IFCs	IFCs for single-cell capture and cDNA preamplification on the C1 system
C1 Single-Cell DNA Seq IFCs	IFCs for single-cell capture and amplification of genomic DNA on the C1 system
C1 Single-Cell mRNA Seq IFCs	IFCs for single-cell capture and generation of cDNA libraries on the C1 system
C1 Single-Cell mRNA Seq HT IFCs	High-throughput (HT) IFCs for single-cell capture and generation of cDNA libraries on the C1 system
Callisto [™] Adherent Cell Culture IFC	IFC for adherent cell culture on the Callisto system
Polaris [™] Single-Cell Dosing mRNA Seq IFC	IFC for single-cell functional genomic studies on the Polaris system, including cell selection, dosing and mRNA sequencing

For technical support visit fluidigm.com/support

 North America
 +1 650 266 6100 | Toll-free (US/CAN): +1 866 358 4354 | techsupport@fluidigm.com
 Latin America
 +1 650 266 6100 | techsupportlatam@fluidigm.com

 Europe/Middle East/Africa/Russia
 +44 1223 598100 | techsupporteurope@fluidigm.com
 China (excluding Hong Kong)
 +86 21 3255 8368 | techsupportchina@fluidigm.com

 Japan
 +81 3 3662 2150 | techsupportjapan@fluidigm.com
 All other Asian countries/India/Australia
 +1 650 266 6100 | techsupportasia@fluidigm.com

For Research Use Only. Not for use in diagnostic procedures.

NOTICE TO PURCHASER: Limited Use License to Perform Preamplification with Fluidigm IFCs

The purchase of this product from Fluidigm Corporation conveys to the purchaser the limited, non-transferable right to the purchaser to use the purchased amount of this product to perform preamplification methods using Fluidigm instruments under license from Thermo Fisher Scientific. No right to resell this product and no other rights (such as real-time PCR methods, apparatus, reagents, or software to perform digital PCR methods) are conveyed by Thermo Fisher Scientific expressly, by implication, or by estoppel. For information on obtaining additional rights, contact outlicensing@lifetech.com or Out Licensing, Thermo Fisher Scientific, 5791 Van Allen Way, Carlsbad, California 92008.

Information in this publication is subject to change without notice. Patent and license information: fluidigm.com/legaInotices. Trademarks: Fluidigm, the Fluidigm logo, C1, Callisto, Open App, and Polaris are trademarks and/or registered trademarks of Fluidigm Corporation in the United States and/or other countries. All other trademarks are the sole property of their respective owners. © 2018 Fluidigm Corporation. All rights reserved. 08/2018