

Juno[™] 96.96 Genotyping Control Line Fluid

Kit Part Number: 100-8574

Intended Use

The purpose of this kit is to provide necessary Juno[™] 96.96 GT Control Line Fluid to support genotyping with TaqMan[®] assays on the Juno[™] system.

Contents

Part Number	Component	Quantity	Storage	Uses
100-8574	Juno [™] 96.96 GT Control Line Fluid	20 syringes	Room temperature	5 IFCs

Handling

Before using reagents, refer to the Genotyping with Juno Getting Started Guide (PN 100-7074).

For technical support visit fluidigm.com/support

EMAIL United States and other countries not in Europe or Asia techsupport@fluidigm.com | Europe techsupporteurope@fluidigm.com | Asia techsupportasia@fluidigm.com
PHONE United States (toll-free) +1 866 358 4354 | Europe +33 1 60 92 42 40 | Japan +81 3 3662 2150 | China +86 21 3255 8368 | All other countries +1 650 266 6100

For Research Use Only. Not for use in diagnostic procedures.

Information in this publication is subject to change without notice. Use standard laboratory safety protocols. Safety data sheet information: fluidigm.com/sds. Patent and license information: fluidigm.com/legalnotices. EU's WEEE directive information: fluidigm.com/compliance. Fluidigm, the Fluidigm logo, and Juno are trademarks or registered trademarks of Fluidigm Corporation in the United States and/or other countries. All other trademarks are the sole property of their respective owners. © 2014 Fluidigm Corporation. All rights reserved. 12/2014