

Maxpar Cell Acquisition Solution

Catalog: 201240

Package size: 200 mL

Storage: Store at 2–8 °C.

Technical Information

Quality Control: Each lot of Maxpar® Cell Acquisition Solution is prepared so as to minimize contamination with heavy elements that could interfere with mass cytometry analysis. Reagent purity is confirmed by elemental analysis using inductively coupled plasma mass spectrometry, verifying that elemental constituents of >100 atomic mass units are below levels that would interfere with CyTOF® mass cytometry assays.

Recommended usage: Maxpar Cell Acquisition Solution is intended for use in preparing cells for acquisition with the WB Injector on the Helios™ system. Refer to the Maxpar Staining with Fresh Fix protocols (PN 400276–400279) for specific instructions regarding usage.

Application notes: Use caution with this reagent to prevent unintentional introduction of elemental contamination, which results in high levels of background signal detection by mass cytometry. Take measures including capping the container when not in use and avoiding contamination from other laboratory reagents.

Description

Maxpar Cell Acquisition Solution is for use in Helios sample acquisition protocols using the WB Injector.

For technical support visit fluidigm.com/support.

North America +1 650 266 6100 | Toll-free (US/CAN): 866 358 4354 | techsupport@fluidigm.com **Latin America** +1 650 266 6100 | techsupportlatam@fluidigm.com
Europe/Middle East/Africa/Russia +44 1223 598100 | techsupporteurope@fluidigm.com **China (excluding Hong Kong)** +86 21 3255 8368 | techsupportchina@fluidigm.com
Japan +81 3 3662 2150 | techsupportjapan@fluidigm.com **All other Asian countries/India/Australia** +1 650 266 6100 | techsupportasia@fluidigm.com

For Research Use Only. Not for use in diagnostic procedures.

Information in this publication is subject to change without notice. **Patent and license information:** fluidigm.com/legalnotices. **Limited Use Label License:** The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased consumable or reagent only with Fluidigm instruments and systems.

Trademarks: Fluidigm, the Fluidigm logo, CyTOF, Helios, and Maxpar are trademarks and/or registered trademarks of Fluidigm Corporation in the United States and/or other countries. © 2018 Fluidigm Corporation. All rights reserved. 12/2018