

Preamp and Reverse Transcription Master Mix— 5 Tubes


Kit Part Number: 100-6301

Purpose



The Preamp and Reverse Transcription Master Mix kit provides necessary reagents to support gene expression using Dynamic Array™ integrated fluidic circuits (IFCs). The kit contains sufficient reagents for 480 reactions.

Contents

PN 100-5581 Preamp Master Mix—5 Tubes

Part Number	Component	Cap Color	Volume/Unit	Unit Quantity	Storage
100-5744	Preamp Master Mix	 Red	106 µL/tube	5 tubes	-20 °C

PN 100-6299 Reverse Transcription Master Mix—5 Tubes

Part Number	Component	Cap Color	Volume/Unit	Unit Quantity	Storage
100-6297	Reverse Transcription Master Mix	 Light green	106 µL/tube	5 tubes	-20 °C
100-5941	PCR Water	 Natural	1.8 mL/tube	1 tube	

Handling

Before using the reagents, see the Real-Time PCR Analysis User Guide (PN 68000088).

Safety

Use standard laboratory safety protocols. Read and understand the safety data sheets (SDSs) before handling chemicals. To obtain SDSs, go to fluidigm.com/sds and search for the SDS using either the product name or the part number.

For technical support visit techsupport.fluidigm.com.

North America +1 650 266 6100 | Toll-free (US/CAN): 866 358 4354 | techsupport@fluidigm.com **Latin America** +1 650 266 6100 | techsupportlatam@fluidigm.com
Europe/Middle East/Africa/Russia +44 1223 598100 | techsupporteurope@fluidigm.com **China (excluding Hong Kong)** +86 21 3255 8368 | techsupportchina@fluidigm.com
Japan +81 3 3662 2150 | techsupportjapan@fluidigm.com **All other Asian countries/India/Australia** +1 650 266 6100 | techsupportasia@fluidigm.com

For Research Use Only. Not for use in diagnostic procedures.

NOTICE TO PURCHASER: Limited Use License to Perform Pre-amplification with Fluidigm IFCs

The purchase of this product from Fluidigm Corporation conveys to the purchaser the limited, non-transferable right to the purchaser to use the purchased amount of this product to perform pre-amplification methods using Fluidigm instruments under license from Thermo Fisher Scientific. No right to resell this product and no other rights (such as real-time PCR methods, apparatus, reagents, or software to perform digital PCR methods) are conveyed by Thermo Fisher Scientific expressly, by implication, or by estoppel. For information on obtaining additional rights, contact outlicensing@lifetech.com or Out Licensing, Thermo Fisher Scientific, 5791 Van Allen Way, Carlsbad, California 92008. Information in this publication is subject to change without notice. **Patent and license information:** fluidigm.com/legalnotices. Fluidigm, the Fluidigm logo, and Dynamic Array are trademarks and/or registered trademarks of Fluidigm Corporation in the United States and/or other countries. © 2018 Fluidigm Corporation. All rights reserved. 12/2018