

# Preamp Master Mix—5 Tubes

**Kit Part Number: 100-5581**

## Purpose

The Preamp Master Mix kit provides a necessary reagent to support gene expression using Dynamic Array™ integrated fluidic circuits (IFCs). The kit contains sufficient reagents for 500 reactions.

## Contents

Part Number	Component	Cap Color	Volume/Unit	Unit Quantity	Storage
100-5744	Preamp Master Mix	 Red	106 µL/tube	5 tubes	-20 °C

## Handling

Before using the reagents, see the Real-Time PCR Analysis User Guide (PN 68000088).

## Safety

Use standard laboratory safety protocols. Read and understand the safety data sheets (SDSs) before handling chemicals. To obtain SDSs, go to [fluidigm.com/sds](http://fluidigm.com/sds) and search for the SDS using either the product name or the part number.

### For technical support visit [techsupport.fluidigm.com](http://techsupport.fluidigm.com).

**North America** +1 650 266 6100 | Toll-free (US/CAN): 866 358 4354 | [techsupport@fluidigm.com](mailto:techsupport@fluidigm.com) **Latin America** +1 650 266 6100 | [techsupportlatam@fluidigm.com](mailto:techsupportlatam@fluidigm.com)  
**Europe/Middle East/Africa/Russia** +44 1223 598100 | [techsupporteurope@fluidigm.com](mailto:techsupporteurope@fluidigm.com) **China (excluding Hong Kong)** +86 21 3255 8368 | [techsupportchina@fluidigm.com](mailto:techsupportchina@fluidigm.com)  
**Japan** +81 3 3662 2150 | [techsupportjapan@fluidigm.com](mailto:techsupportjapan@fluidigm.com) **All other Asian countries/India/Australia** +1 650 266 6100 | [techsupportasia@fluidigm.com](mailto:techsupportasia@fluidigm.com)

### For Research Use Only. Not for use in diagnostic procedures.

**NOTICE TO PURCHASER:** Limited Use License to Perform Pre-amplification with Fluidigm IFCs

The purchase of this product from Fluidigm Corporation conveys to the purchaser the limited, non-transferable right to the purchaser to use the purchased amount of this product to perform pre-amplification methods using Fluidigm instruments under license from Thermo Fisher Scientific. No right to resell this product and no other rights (such as real-time PCR methods, apparatus, reagents, or software to perform digital PCR methods) are conveyed by Thermo Fisher Scientific expressly, by implication, or by estoppel. For information on obtaining additional rights, contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Thermo Fisher Scientific, 5791 Van Allen Way, Carlsbad, California 92008. Information in this publication is subject to change without notice. **Patent and license information:** [fluidigm.com/legalnotices](http://fluidigm.com/legalnotices). Fluidigm, the Fluidigm logo, and Dynamic Array are trademarks and/or registered trademarks of Fluidigm Corporation in the United States and/or other countries. © 2018 Fluidigm Corporation. All rights reserved. 12/2018