

## Skimmer Reducer Assembly

Catalog#: 101802  
Package Size: 1 unit  
Instrument: CyTOF and CyTOF 2

### Technical Information

**Description:**

The CyTOF® Skimmer-Reducer assembly comprises the last two parts of the 3-cone interface which transmits ionized clouds generated in the plasma to the low pressure mass analyzer chambers of the cytometer. The Skimmer and Reducer are made of nickel and have a 1mm and 1.2mm diameter aperture, respectively. The assembly is installed behind the Sampler Cone.



### Contact Information:

Sales: [sales@DVSSciences.com](mailto:sales@DVSSciences.com) | Support: [Support.NorthAmerica@DVSSciences.com](mailto:Support.NorthAmerica@DVSSciences.com) or [Support.Europe@DVSSciences.com](mailto:Support.Europe@DVSSciences.com)  
[www.DVSSciences.com](http://www.DVSSciences.com) | For assistance by phone: +1-855-DVS-CYTO

**This product is intended for Research Use Only. This product is not intended for therapeutic or diagnostic purposes in humans or animals.**

**Vendor info & description quote as needed.**

List supplier of raw material antibody here only if targets/clones are unique, licensed, and/or indicated as required in our agreement.

© 2011 DVS Sciences Inc. The trademarks mentioned herein are the property of DVS Sciences Inc. or their respective owners.

DVS Sciences Inc. products are warranted to meet stated product specifications and to conform to label descriptions when used and stored properly. Unless otherwise stated, this warranty is limited to 6 months from date of sales for products used, handled and stored according to DVS Sciences, Inc. instructions. DVS Sciences, Inc. sole liability for the product is limited to replacement of the product or refund of the purchase price. Product may not be resold, modified for resale or used to manufacture commercial products without prior written approval from DVS Sciences Inc.

