


20X GE Sample Loading Reagent—5 Tubes

Kit Part Number: 100-7610

Purpose

The 20X GE Sample Loading Reagent kit provides a necessary reagent to support gene expression using Dynamic Array™ IFCs. The kit contains sufficient reagents for fifty 48.48 integrated fluidic circuits (IFCs) or twenty-five 96.96 IFCs.

Contents

Part Number	Component	Cap Color	Volume/Unit	Unit Quantity	Storage
85000735	20X GE Sample Loading Reagent	 Blue	250 µL/tube	5 tubes	-20 °C

Handling

Before using this kit, see the Real-Time PCR Analysis User Guide (PN 68000088).

Safety

Use standard laboratory safety protocols. Read and understand the safety data sheets (SDSs) before handling chemicals. To obtain SDSs, go to fluidigm.com/sds and search for the SDS using either the product name or the part number.

For technical support visit techsupport.fluidigm.com.

North America +1 650 266 6100 | Toll-free (US/CAN): 866 358 4354 | techsupport@fluidigm.com **Latin America** +1 650 266 6100 | techsupportlatam@fluidigm.com
Europe/Middle East/Africa/Russia +44 1223 598100 | techsupporteurope@fluidigm.com **China (excluding Hong Kong)** +86 21 3255 8368 | techsupportchina@fluidigm.com
Japan +81 3 3662 2150 | techsupportjapan@fluidigm.com **All other Asian countries/India/Australia** +1 650 266 6100 | techsupportasia@fluidigm.com

For Research Use Only. Not for use in diagnostic procedures.

Information in this publication is subject to change without notice. **Patent and license information:** fluidigm.com/legalnotices. Fluidigm, the Fluidigm logo, and Dynamic Array are trademarks and/or registered trademarks of Fluidigm Corporation in the United States and/or other countries. © 2018 Fluidigm Corporation. All rights reserved.
12/2018