


# 2X Assay Loading Reagent—5 Tubes

**Kit Part Number: 100-7611**

## Purpose

The 2X Assay Loading Reagent kit provides a necessary reagent to support gene expression and genotyping using Dynamic Array™ IFCs. The kit contains sufficient reagents for fifty 48.48 integrated fluidic circuits (IFCs) or twenty-five 96.96 IFCs.

## Contents

Part Number	Component	Cap Color	Volume/Unit	Unit Quantity	Storage
85000736	2X Assay Loading Reagent	 Yellow	1.5 mL/tube	5 tubes	-20 °C

## Handling

Before using this kit, see the Real-Time PCR Analysis User Guide (PN 68000088) or the SNP Genotyping Analysis User Guide (PN 68000098).

## Safety

Use standard laboratory safety protocols. Read and understand the safety data sheets (SDSs) before handling chemicals. To obtain SDSs, go to [fluidigm.com/sds](http://fluidigm.com/sds) and search for the SDS using either the product name or the part number.

### For technical support visit [techsupport.fluidigm.com](http://techsupport.fluidigm.com).

**North America** +1 650 266 6100 | Toll-free (US/CAN): 866 358 4354 | [techsupport@fluidigm.com](mailto:techsupport@fluidigm.com) **Latin America** +1 650 266 6100 | [techsupportlatam@fluidigm.com](mailto:techsupportlatam@fluidigm.com)  
**Europe/Middle East/Africa/Russia** +44 1223 598100 | [techsupporteurope@fluidigm.com](mailto:techsupporteurope@fluidigm.com) **China (excluding Hong Kong)** +86 21 3255 8368 | [techsupportchina@fluidigm.com](mailto:techsupportchina@fluidigm.com)  
**Japan** +81 3 3662 2150 | [techsupportjapan@fluidigm.com](mailto:techsupportjapan@fluidigm.com) **All other Asian countries/India/Australia** +1 650 266 6100 | [techsupportasia@fluidigm.com](mailto:techsupportasia@fluidigm.com)

### For Research Use Only. Not for use in diagnostic procedures.

Information in this publication is subject to change without notice. **Patent and license information:** [fluidigm.com/legalnotices](http://fluidigm.com/legalnotices). Fluidigm, the Fluidigm logo, and Dynamic Array are trademarks and/or registered trademarks of Fluidigm Corporation in the United States and/or other countries. © 2018 Fluidigm Corporation. All rights reserved.  
12/2018